

# Size calculation

## Elbow Orthosis Grade 2



### Measurement

During measuring the animal is


- ☐ Standing / weight bearing - recommended!
- ☐ Lying down

### Circumference

1. Measure the circumference of the paw \_\_\_\_\_
2. Measure the circumference just above the paw \_\_\_\_\_
3. Measure the circumference at the carpal joint \_\_\_\_\_
4. Measure the circumference midway between the elbow and carpal \_\_\_\_\_
5. Measure the circumference at elbow \_\_\_\_\_
6. Measure the circumference at the maximum proximal region/ armpit \_\_\_\_\_
7. Measure from the center of the line from the armpit across to the same point on the opposite side.  
\_\_\_\_\_

### Length

1. Measure from the carpal joint to the foot pad \_\_\_\_\_
2. Measure from the carpal joint to the elbow \_\_\_\_\_
3. Measure from elbow to the measurement on # 6  
\_\_\_\_\_



Send measurements to Scandi Orthopedic.